Lock

Bipolar head
Bipolarna glava

**PRODUCT DESCRIPTION
AND SURGICAL TECHNIQUE
OPIS PROIZVODA
I KIRURŠKA TEHNIKA**

Index

Sadržaj

pag.	3	HISTORY AND PHILOSOPHY	<i>POVIJEST I FILOZOFIJA</i>
pag.	4	MATERIALS	<i>MATERIJALI</i>
pag.	5	BIOMECHANICAL CHARACTERISTICS	<i>BIOMEHANIČKE KARATERSTIKE</i>
pag.	6	PRODUCT CHARACTERISTICS	<i>KARAKTERISTIKE PROIZVODA</i>
pag.	7	INDICATIONS	<i>INDIKACIJE</i>
pag.	7	SIZES	<i>VELIČINE</i>
pag.	8	PRODUCT CODES	<i>ŠIFRE PROIZVODA</i>
pag.	10	INSTRUMENT SET	<i>SET INSTRUMENATA</i>
pag.	12	'LOCK BIPOLAR HEAD' SURGICAL TECHNIQUE	<i>KIRURŠKA TEHNIKA</i> <i>"LOCK BIPOLARNE GLAVE"</i>
pag.	12	1. Sizing the bipolar head	<i>1. Određivanje veličine bipolarne glave</i>
pag.	13	2. Trial reduction	<i>2. Probna redukcija</i>
pag.	14	3. Bipolar head assembling	<i>3. Montaža bipolarne glave</i>
pag.	16	4. Final reduction and removal of the components	<i>4. Konačna redukcija i uklanjanje drugih komponenti</i>
pag.	17	BIBLIOGRAPHY	<i>BIBLIOGRAFIJA</i>
pag.	17	ARTICLES	<i>ČLANCI</i>

History and Philosophy

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The Lock bipolar head proposes the philosophy of the partial hip prosthesis with forefront criteria of design, materials and ergonomics.

An highly congruent metallic polished shell working on the centers self-aligning principle, is equipped by conical UHMWPE liner and a simple retentive ring mechanism.

An essential instrumentation set allows to perform a fast and effective surgical technique.

Lock bipolarna glava pripada filozofiji parcijalne proteze kuka s naprednim kriterijem dizajna, materijala i ergonomije.

Vro kongruentna metalna polirana školjka koja radi na principu automatskog poravnavanja, opremljena je koničnim UHMWPE umetkom i jednostavnim prstenastim mehanizmom.

Osnovni set instrumenata omogućava brzu i učinkovitu kiruršku tehniku.

Materials

Materijali

The Lock bipolar heads are available in AISI 316L stainless steel or CoCrMo alloy, well known materials for the tribological characteristic.

The polished surface finish and the high congruence to the acetabular cavity allows a natural fluency and adaptation to the chondral surface.


Each shell is uniquely combined with a liner and a retentive ring all made of UHMWPE.

Sterilization of the items is provided by means of ethylene oxide (Eto).

This method preserves intact the characteristics of the UHMWPE, because it only acts at the surface level of the solid parts, without long term degrading effects for the tribological and mechanical properties.

Physical chemical interactions are not possible by the sterilizing agent, therefore the polymer cannot be affected by structural alterations as per irradiating processes.

Bipolarne glave s vratom na raspolaganju su u izvedbi iz AISI 316L nehrđajućeg čelika ili CoCrMo legure, dobro poznati materijali za tribološke karakteristike.

Polirana površina i visoka kongruencija acetabularne šupljine, omogućava prirodni tok i adaptaciju hondralne površine.

Svaka školjka je jedinstveno kombinirana s umetkom i retentivnim prstenom, sve u UHMWPE izvedbi.

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
A nearly hemispheric profile characterizes the outer shell. It increases the R.O.M with respect to the inner cephalic modular component, and decreases the risk ratio of accidental peripheral contacts with the femoral neck.

Skoro hemisferični profil karakterizira vanjsku školjku. Time se povećava R.O.M u odnosu na unutarnju cefaličku modularnu komponentu, te ujedno i smanjuje postotak rizika slučajnog perifernog kontakta s femoralnim vratom.

The positive eccentricity of the inner articulating head allows the self-alignment of the two centers of rotation with respect to the bearing load direction in every mechanical condition.

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
A distributed progression for the shell sizes, facilitates the anatomical choice.


Distribuirana progresija veličine školjke olakšava anatomske izbor.

Moreover, 6 different lengths of the inner head allow the user to re-establish the optimal tissue balancing.

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



Conical coupling.
Konični spojevi

Hemispherical coupling
Hemisferni spojevi

The inner head size is of 28mm of diameter.

Such a geometry allows a range of motion beyond 75° with respect to the equatorial edge of the outer shell (NB: possible values with Lima-Lto femoral stems).

The mechanical connection between liner and shell is self-locking and possible by Morse cone design. This method, originally introduced in 1978 by Lima-Lto, totally avoids any risks of micro-movement and consequent back side wear phenomena at the liner-shell interface, reported for the snap hemispherical coupling.

Shell, liner and ring are in unique packaging and suitable for a fast intra-operative assemblage. Sizing and implanting can be performed with the support of a simple and efficient instrumentation.

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Takva geometrija omogućava domet kretanja veći od 75° u odbosu na ekvatorni rub vanjske školjke (NB: moguće vrijednosti s Lima-Lto femoralnim stemom).

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Indications

Traumatic pathology of the head and the femoral neck in elder patients.

Degenerating pathologies of the hip with peculiar effect mostly to the femoral head such as a-vascular necroses disease.

Special use for selected case of revision hip arthroplasty whereas a serious acetabular deficiency still allow the maintenance of an adequate articulation congruency.

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Sizes

The Lock bipolar heads are available in 14 sizes, increasing of millimeter in millimeter, from the diameter 41mm to the diameter 54mm. Upon request 3 bigger sizes of diameters 55mm, 56mm and 57mm.

Those shell are combined with 28mm inner articulation.

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
'Lock Bipolar Head' Surgical Technique

Lock bipolarna glava - Kirurška tehnika

[1]

Using the desired surgical technique, expose the acetabulum so as to view it adequately, then remove possible osteophytes and soft tissues surrounding the cotyloid rim.

1. Sizing the bipolar head

A suitable size for the bipolar head can be chosen by comparing the main overalls of the articular ends.

[2]

In fact, if the head of the femur is still in good conditions, a femoral caliper can be used to size its diameter (Fig. 1-2).

A trial bipolar head (Fig. 3-4) can also be used by a special inserter to size the acetabular site (Fig. 5-6).

[3]

To make it a fast maneuver of sizing the inserter is designed to offer two clamping ends.

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[5]

2. Trial reduction

Insert a trial femoral head onto the femoral stem (use only Lima Lto femoral prosthetic implantable devices) then insert the trial bipolar head provide to reduce the articulation (Fig. 7-8).

[6]

The trial bipolar heads are designed to self-retain the femoral trial heads, as enough as to allow a feasible and reliable trial reduction, without altering the range of motion and tissue balancing.

The joint is reduced.

[7]

Check the length of the trial femoral head is adequate to compensate tissue balancing and range of motion in the way to avoid the dislocation (Fig. 9).

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3. Bipolar head assembling

Read on the packaging label the sizes of the chosen component.

Metal shell and liner. Preserve apart the retentive ring.

Insert the liner by hand into the metal shell (Fig. 10)

Using a gauze during the maneuver (Fig. 11) adjust the liner in deepest position with the pusher liner by a slight pressure (Fig. 12).

Warning: The first maneuver is always insert the inner femoral head onto the femoral stem, than hang the metal shell onto the femoral head with the snapping ring (use only Lima Lto femoral prosthetic implantable devices).

Having already coupled the inner femoral head to the metal shell it is very dangerous to hang the bipolar head onto the femoral stem taper. In fact there is a risk of misalignment because it is not possible to keep the orientation of inner head taper relatively to the direction of maneuver

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The ring should already have been connected on the femoral neck before the femoral head (Fig. 13). Check that the two holes at the chink extremity of the ring are towards the femoral side and not towards the acetabulum (Fig. 14).

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If the head were ceramic (mind to clean and dry the neck taper) fix it by pushing and rotating it along the neck axis (Fig. 15-16).

Ukoliko su glave bile keramičke (pazite da je vratni konus očišćen i osušen) fiksirajte ih tako da ih gurnete i rotirate uzduž vratne osi. (sl. 15-16).

The metal shell is coupled onto the head and the ring is drawn. Push the shell against the head and at the same time pull the ring until obtaining a fast elastic expansion, which indicates the ring snapped into the bipolar seat (Fig. 17).

Metalna školjka spojena je na glavu a prsten je povučen. Gurnite školjku prema glavi i istodobno vucite prsten sve dok niste doбили elastičnu ekspanziju koja ukazuje na to da je prsten skočio u bipolarno sjedište (sl. 17).

[18]

4. Final reduction and removal of the components

Reduce the articulation and check for the correct ratios (Fig. 18).

If necessary, the various prosthetic components can be disassembled and removed. The metallic shell can be disconnected from the inner femoral head by removing the retentive ring with the pliers (Fig. 19-20).

To remove the inner femoral head, just beat the head base axially with a punch (not provided with the instrument set).

Warning!

In case of head removal only, never use a ceramic head on the same taper. In this case, the stem must be replaced as well.

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Articles

Članci

A Bipolar Cup with Ceramic Insert for the Hip Degenerative Pathology: Early Experiences

Bipolarna čašica s keramičkim umetkom za degenerativne patologije kuka: rana iskustva

L.Fantasia - D.Cornacchia -

M. Agamennone - P.La Floresta

2nd BIOLOX® Symposium, Proceedings 1997

UHMWPE for Arthroplasty

UHMWPE za artroplastiku

L. Costa - E.M. Brach del Prever

Ed. Minerva Medica

[illegible]

For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

Glede daljnjih informacija o našim proizvodima molimo da pogledate našu web-stranicu www.lima.it.

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Lima Deutschland GmbH