Logica

Femoral Stems
Steli femorali

Standard

Mirror

CS

**SURGICAL TECHNIQUE
TECNICA OPERATORIA**

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
Product Codes

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Cemented STANDARD Stem
Stelo cementato STANDARD

Cemented MIRROR Stem
Stelo cementato MIRROR

Cemented MIRROR Stem
Stelo cementato MIRROR

Cemented STANDARD Stems with Centralizer - Taper 12/14 Steli cementati STANDARD con centratore - Cono 12/14

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
1. Pre-operative Planning

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Logica Mirror

In order to obtain the best results it is highly recommended to make a pre-operative planning by using X-rays, which give a 15% enlargement of profiles. Use high quality frontal X-rays with a proper contrast, including the whole length of stems on pre-operative X-rays.

Determine the stem size and the neck resection level so that the great trochanter upper level coincides with the medium femoral head centre (Fig. 1-2).

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
2. Position of the Patient

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[3]

Note. The first figure of surgical technique shows the posterolateral approach with the patient in lateral decubitus position. The next technique steps also refer to the supine position and any surgical access.

Patient in lateral decubitus position. Posterolateral incision (Fig. 3). With the fascia lata opened and the lateral rotators sectioned, the articular capsule is incised and then the femoral epiphysis is dislocated. The femoral head is exposed.

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3. Neck Resection

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[4]

The transcervical neck resection is carried out as planned by means of an oscillating saw, with an angle of 45° with respect to the diaphyseal axis (Fig. 4).

Then proceed to the acetabular bone preparation and corresponding replacement.

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4. Preparation of the Medullary Canal

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Canal Chisel
Fustella per canale diafisario
cod. 9095.10.160

[5]

The diaphysis is perforated along the femoral axis by means of the **canal chisel (G)** (Fig. 5-6).

[6]

The instrument set includes two different rasp handles. If the **rasp handle (B)** - suitable for posterior approach and with the patient in lateral decubitus position - is used, open the lever and insert the rasp bolt through the handle, with the medial side oriented towards the lever. Then tighten the lever until it is closed (Fig. 7-8).

The canal is then prepared with the **rasps (A)**, starting from the smallest rasp and maintaining the correct antiversion (approx. 15°).

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

Rasp
Raspa
cod. 9037.30.002...

[9]


The **clamp rasp handle (C)** can be used likewise, mainly with obese patients or for direct lateral approaches and with the patient in the supine position.

With the handle lever engage the proximal tang of the rasp at the handle end. Once the rasp is coupled, push the locking lever against the handle (Fig. 10).


Note. The clamp teeth can be separated to wash it thoroughly (Fig. 11-12).

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[12]

[13]

Finish the medullary canal preparation by filling the epiphysis seat properly. Stop the rasp insertion next to the cephalic resection or one of the three points carved on the rasp, based on the insertion level chosen for the prosthesis (Fig. 14).

Open the closing lever and remove the rasp handle, leaving the rasp inside the medullary canal.

Note. The same rasp series can be used for all Logica system stems.

[14]

When inserting Standard or Mirror cemented stems, the seat created keeps the stem in a stable position enclosed by a level cement layer with a minimum thickness of 1.8 mm. If the cement layer must be reduced, the final stem with one size higher than the last rasp used can be inserted.

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Trial Head
Testa di prova
cod. 9095.10.340...

[15]

Keeping the rasp in situ, it is now possible to determine the size of the femoral head to be used.

Insert the **trial femoral heads (D-I)** (Fig. 15-16) into the rasp neck and reduce the articulation by housing the trial femoral head in the acetabular cup already inserted (Fig. 17).

Select the length of the femoral head to be used (S, M, L) based on the articular tension.

Then remove the trial head and the rasp by means of the rasp handle.

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[18]

FOR STANDARD AND MIRROR STEMS:

After installing a distal plug (available on request), fill the prosthesis seat with acrylic cement backwards (Fig. 18).

Note. LIMA-LTO supplies cementation and pressurisation kits together with femoral tampons, for an optimal and efficient cementation technique.

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
[19]

FOR STANDARD STEMS:

Insert the stem manually until the collar seat is in contact with the resection level. During this operation, the collar facilitates the cement pressurisation on the proximal area.

Keep the stem in a stable position until the polymerisation process is terminated by means of the **impactor stem positioner (F)**.

Remove the cement excess leakage before the polymerisation process is finished.

FOR MIRROR STEMS:

Introduce the stem manually up to the required depth indicated by the three different insertion levels (Fig. 20). Keep the stem in a stable position until the polymerisation process is terminated by means of the **impactor stem positioner (F)** (Fig. 21).

Remove the cement excess leakage before the polymerisation process is finished.

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
[22]

FOR CS STEMS:

Screw the **impactor-extractor (E)** into the threaded hole on the stem shoulder, then proceed to its introduction by progressive percussions.

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[25]

8. Final Trial Reduction

Insert the trial femoral heads (available on request - Fig. 25) through the stem taper and reduce the articulation by housing the head in the acetabular cup previously inserted.

Select the length of the femoral head to be used (S, M, L) based on the articular tension (Fig. 26).

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9. Head Insertion

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[27]

In the sterile box take the head with the diameter and length previously determined.

Before the head insertion, clean and dry the taper neck carefully. This operation is essential, especially when using ceramic heads.

Secure the head in the neck by pushing it along the neck axis and rotating it simultaneously (Fig. 27).

Remove the protective cap (Fig. 28).

[28]

Then hit slightly with a plastic beater (available on request) along the taper axis (Fig. 29-30).

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10. Final Reduction

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[31]

Then reduce the articulation after having checked that the articular surfaces are perfectly clean (Fig. 31).

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11. Removal of the Components

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
If needed, all components can be removed. To remove the femoral head only, just hit the head base axially with a punch.

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
Warning! When replacing a ceramic head, never use a second ceramic head on the same stem. In this case, the stem must be replaced as well.

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The stem can be removed with the **impactor-extractor (E)** (Fig. 32) by screwing it into the threaded hole on the stem shoulder (Fig. 33-34).

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[illegible]

[illegible]

For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

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